P0#2B1058 LS802-B 说明书 (A0)

印色: 单黑 + Pantone 186 c

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(P0#2B0846 改:P4/P37)

TENSIOSCREEN

USER MANUAL | GUIDE D'UTILISATION

NL Handleiding | **IT** Manuale di istruzioni | **ES** Manual de instrucciones

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www.terrailon.com

Made for / Compatibilité


iPhone® 4S/5/5C/5S/6/6Plus
iPod® Touch 5th generation

iPad® 3/4/Air
iPad® Mini

Android 4.3
Bluetooth Smart 4.0 / Bluetooth Smart Ready

CE 0123

iPhone is registered trademark of Apple Inc.

Android is a trademark of Google Inc.

The Bluetooth® word mark and logos are registered trademarks
owned by Bluetooth SIG, Inc.

Other trademarks and trade names are those of their respective owners.

Tensioscreen

EN	TERRAILLON CONNECTED BLOOD PRESSURE MONITOR INSTRUCTION MANUAL	4
FR	TENSIOMÈTRE CONNECTÉE TERRAILLON NOTICE D'UTILISATION	8
NL	TERRAILLON BLOEDDRUKMETER HANDLEIDING	12
IT	SFIGMOMANOMETRO TERRAILLON MANUALE DI ISTRUZIONI	16
ES	TENSIÓMETRO TERRAILLON MANUAL DE INSTRUCCIONES	20
DE	BLUTDRUCKMESSGERÄT TERRAILLON BEDIENUNGSANLEITUNG	24
PT	MONITOR DE PRESSÃO SANGUINEA TERRAILLON MANUAL DE INSTRUÇÕES	28

Thank you very much for selecting the Terraillon Arm Blood Pressure Monitor TENSIOSCREEN.



This unit allows you to monitor your blood pressure. This device is not intended to be a diagnostic device. Always contact your doctor for advice, diagnostic and treatment.

FIND MORE INFORMATION ON

<http://www.terraillon.com>

PRODUCT OVERVIEW


- SYS** Systolic Blood Pressure
- DIA** Diastolic Blood Pressure
- Pul/min** Pulse
- mmHg** Unit
- (☎)** Data transmitting
- (🔋)** Low Battery
- (📄)** Data pending to transmit
- ⌚** Time [Hour/Minute - Month/Day]
- ❤️** Irregular Heartbeat Detector

FIRST USE

- A** Insert the batteries (4 x AA) and press and hold User1 button to set time. Press button **Q** to change the numeral then press User1 button to confirm. After confirming [TIME] and [DATE], the LCD will display "Done".
- B** Download the Terraillon Wellness Coach application.


- C** Create your account.


- D** Turn on Bluetooth on your Smartphone (Setting > Bluetooth > ON).
- E** Press and hold User2 button to start pair-up.

F

START MEASURING

- A** Tie the cuff
- B** Press User1 or User2 button to turn on the monitor, and it will complete the measurement process automatically

RECORDS

Press button  to access the memory. Press User1 or User2 to recall the measurement history of each user.

Press and hold button  for 3 seconds to clear the memory. When the LCD display "dEL ALL", press button  to confirm.

TIPS FOR MEASUREMENT

	Optimal	Normal	High-normal	Mild	Moderate	Severe
SYS	<120	120-129	130-139	140-159	160-179	≥180
DIA	<80	80-84	85-89	90-99	100-109	≥110

* The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH).

TROUBLESHOOTING

- E2** The cuff is very tight.
- E3** The pressure of the cuff is excess.
- E9** Product has not been activated.
- E10/11** The monitor detected motion while measuring.
- E20** Measurement process does not detect the pulse signal.
- E21** Measure incorrectly.
- Eexx** A calibration error occurred.
- Lo** Low Battery. Replace all four batteries in the same time.

WARRANTY & ENVIRONMENT PROTECTION

This product is guaranteed for 2 years against material and manufacturing defects. During this period, such defects will be repaired free of charge (proof of purchase must be presented if the scale is under warranty). This warranty does not cover damage resulting from accidents, improper use or negligence. If you have a complaint please contact the store where you purchased the product first.

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

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Safety Information

The below signs might be in the user manual, labeling or other components. They are the requirement of standard and using.

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol for "TYPE BF APPLIED PARTS"
CE 0123	Symbol for "COMPLIES WITH MDD93/42/EEC REQUIREMENTS"		Symbol for "ENVIRONMENT PROTECTION – Waste electrical products should not be disposed of with household waste. Please follow local guidelines"
	Symbol for "MANUFACTURER"		
SN	Symbol for "SERIAL NUMBER"		Symbol for "Authorised Representative in the European Community"
	For indoor use only		Symbol for "Class II Equipment"
	The Bluetooth Combination Mark	F1	T1A/250V Ø3.6*10CCC


CAUTION

- * Please do read this user manual carefully and thoroughly before use.
- * This device is intended for adult use only.
- * This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the upper arm or for functions other than obtaining a blood pressure measurement.
- * Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Please start or end medical treatment basing solely on physician's treatment advice.
- * If you are taking medication, consult your physician to determine the most appropriate time for your measurement. Never change a prescribed medication without your physician's consent.
- * This unit is not suitable for continuous monitoring during medical emergencies or operations.
- * If the pressure of the cuff exceeds 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when its pressure exceeds 40 kPa (300 mmHg), detach the cuff from the upper arm and press the homologous button to stop inflation.
- * Please do read this user manual carefully and thoroughly before use.
- * Do not use the monitor under the conditions of strong electromagnetic field (e.g. mobile) that radiates interference signal or electrical fast transient / burst signal, especially when the AC adaptor is applied.
- * Do not touch the output of AC adaptor and the patient simultaneously.
- * The device is not AP/APG equipment. It is not suitable for use in the presence of a flammable anesthetic mixture with air (or oxygen, nitrous oxide).
- * Please keep the unit out of reach of infants or children, since inhalation or swallowing of small parts is dangerous or even fatal.
- * Please use ACCESSORIES and detachable parts specified / authorised by MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user / patient.
- * The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5: 2009 and ISO 10993-10:2010. It will not cause any potential allergic reaction or contact injury.
- * Please make sure the unit functions safely and it is in proper working conditions before use.

Tips for Measurement

It can cause inaccuracy if the measurement is taken in the following circumstances.

Within 1 hour
after dinner or drinking

Immediate measurement
after tea, coffe, smoking

Within 20 minutes
after taking a bath

When talking or moving
your fingers

In a very cold
environemnt


When you want
to discharge urine


Maintenance


To obtain the best performance, please follow below instructions.


Put in a dry place
and avoid the sunshine


Avoid immersing it in the water.
Clean it with a dry cloth in case


Avoid shaking and collision

Avoid dusty environment and
unstable temperature surrounding

Use the slightly damp cloth
to remove the dirt

Avoid washing the cuff

Cleaning: Before Use - Pick out the whole unit of the storage bag. Use the soft cloth to remove the dirt on the monitor and apply some alcohol to disinfect the cuff before tying the cuff.

After Use - Use the soft cloth to wipe the unit and apply some alcohol to disinfect the cuff before putting the whole unit back in the bag.


Please always disinfect the cuff before applying to another patient.

Please follow the instructions for correct replacement of interchangeable or detachable parts specified by SERVICE PERSONNEL of MANUFACTURER as "replaceable".

Disposal: Degraded sensors may result in inaccurate measurement while loosened electrodes may cause the monitor's failure to power on. The expected life of the monitor is two years. Please dispose of ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.


What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, which is called systolic pressure. When the ventricles relax, the blood pressure reaches its minimum value in the cycle, which is called diastolic pressure.

What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

CAUTION

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

Level	Óptima	Normal	Normal-alta	Leve	Moderada	Grave
Blood Pressure (mm Hg)						
SYS	<120	120-129	130-139	140-159	160-179	≥180
DIA	<80	80-84	85-89	90-99	100-109	≥110

Irregular Heartbeat Detector

This Blood Pressure Monitor is equipped with an intelligent function of Irregular Heartbeat (IHB) Detector. During each measurement, this equipment records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, this equipment will light up the IHB symbol on the screen when displaying the measuring result.


CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Why my blood pressure is varies even in one day?

1. Individual blood pressure varies every in one day, it also affected by the way you tie your cuff and the your measurement position, so please take the measurement at the same condition.
2. The varies of the pressure is greater if the person take medicine.
3. Waiting at least 4-5 minutes for another measurement.

Why the blood pressure I get from the hospital is different from home?


The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the "white coat" in hospital which makes the results are higher than the ones at home.

The result is different that you measuring on the different arm.

It is ok for both arms, but there will be some different results for different arm, so suggest you measure the same arm every time.

The attention need to pay when you measure you blood pressure at home:

- If the cuff is tied properly.
- If the cuff is too tight or too loose.
- If the cuff is tied on the wrist.
- If you feel anxious pressured.
- You had better take deep breath 2-3 times before beginning.
- Advice: adjust yourself for 4-5 minutes until you calm down.

SPECIFICATIONS

Power supply	Battery powered mode: 6VDC 4*AA-size alkaline batteries AC adaptor powered mode: 100-240 V~, 50-60 Hz, 400 mA (Can be supplied by AC adaptor model UE08WCP-060100SPA only!)
Display mode	Digital LCD V.A. = 78mm x 92mm
Measurement mode	Oscillographic testing modev
Measurement range	Pressure: 0-40kpa(0~300mmHg) pulse value: (40-199)times/minute
Accuracy	Pressure: 5°C-40°C withi n±0.4kpa(3mmHg) pulse value: ±5%
Working condition	Temperature:5°C-40°C Relative Humidity ≤85%RH Atmospheric Pressure: 86-106 kPa
Storage & transportation condition	Temperature:-20°C-60°C Relative Humidity 10%-93%RH Atmospheric Pressure: 50-106 kPa
Measurement perimeter of the upper arm	About 22cm-42cm
Weight	Approx.385g(Excluding the dry cells)
External dimensions	Approx.120*160*69mm
Attachment	4*AA alkaline batteries,user manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor Powered Mode: Class II ME Equipment
IP Classification	IP22
Software Version	V01

WARNING: No modification of this equipment is allowed.

Complied European Standards List

Risk Management	EN/ISO 14971:2007
Labeling	EN 980:2008
User Manual	EN 1041:2008
General Requirements for Safety	EN 60601-1:2006/A1:2012 EN 62304:2006/AC:2008 EN 60601-1-6:2010 EN 60601-1-11:2010
Non-invasive Sphygmomanometers General Requirements	EN 1060-1:1995+A2:2009 EN 1060-3:1997+A2:2009 EN 1060-4:2004
Electromagnetic Compatibility	EN 60601-1-2:2007/AC:2010

EMC Guidance

Table 1 Guidance and manufacturer's declaration – electromagnetic emissions for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission		
The LS802-B is intended for use in the electromagnetic environment specified below. The customer of the user of the LS802-B should assure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions C ISPR 11	Group 1	The LS802-B uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	The LS802-B is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

**Table 2 – Guidance and MANUFACTURER’S declaration – electromagnetic IMMUNITY –
for all ME EQUIPMENT and ME SYSTEMS**

Guidance and manufacturer’s declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Table 4 – Guidance and MANUFACTURER’S declaration – electromagnetic IMMUNITY – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer’s declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.167 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1.167 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.333 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.</p> <p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 6 – Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.167 \sqrt{P}$	80 MHz to 800 MHz $d = 1.167 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.333 \sqrt{P}$
0,01	0.117	0.117	0.233
0,1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.690	3.690	7.378
100	11.67	11.67	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance <i>d</i> in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

GUANGDONG TRANSTEK MEDICAL ELECTRONICS CO., LTD
 Zone A, No.105 ,Dongli Road, Torch Development District,
 Zhongshan,528437,Guangdong,China

MDSS - Medical Device Safety Service GmbH
 Schiffgraben 41, 30175 Hannover, Germany

Terraillon SAS France & Headquarters

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Terraillon Asia Pacific Ltd

4/F, Eastern Centre
1065 King's Road
Quarry Bay - Hong Kong
Tel: + 852 (0)2960 7200
customerservice@terraillon-asia.com

Terraillon Corp USA

contact@terraillon.fr

Terraillon®

www.terraillon.com

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Terraillon UK Ltd