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EN

ULTRAFAST THERMOMETER

Before use, read the instructions and warnings given in this handbook carefully and keep the handbook in a safe place. Digital thermometer with a large LCD display for easy reading. Fast and accurate fixed with a sound signal when the temperature reading has been taken. Stores the last read temperature, to evaluate any variations. Extremely easy to use, it is kept in a practical case to avoid risk of breaking. The flexible probes can be washed, and it turns off automatically.

GENERAL WARNINGS
1- Self-measuring means control, not diagnosis and treatment, any unusual readings must be discussed with your doctor.
2- Under no circumstances must you charge the medicine doses that have been prescribed by your doctor.
3- Do not squeeze, bend, drop or break the thermometer.
4- Keep it out of the reach of children.
5- Do not expose to sunrays or direct sources of heat. Keep in a dry clean place.
6- Do not immerse the thermometer in water or other liquids, only the probe is waterproof.
7- This thermometer reads only the armpit, rectal or oral temperature. Do not use it for reading armpit temperature.
8- In certain environments subject to radio interference, the product may operate incorrectly, in such cases, move the thermometer away from this area, or remove the battery and replace it.

INSTRUCTIONS FOR USE
1- Press the ON/OFF key. For a few seconds all the function symbols light up on the display. The thermometer will then display 36.5°C then L, and the °C symbol will flash. The thermometer is now ready to use. 2- Place the probe as follows: IN THE MOUTH: place the probe under the tongue and keep the mouth closed for about 10-15 seconds, until the temperature stabilises. The standard temperature for this type of reading is 37°C UNDER THE ARM: place the probe under the armpit and keep the arm well pressed against the chest for about 10-15 seconds, until the temperature stabilises. The standard temperature for this type of reading is 36.7°C IN THE RECTUM: lubricate the probe with a water-soluble gelatinous solution. Do not use Vaseline. Insert the probe into the rectum for about 1 cm., until the temperature stabilises (about 10-15 seconds). The standard temperature for this type of reading is 37.0°C.
3- When the thermometer gives off the sound signal, remove it. The reading is terminated and shown on the display. The thermometer turns off automatically after 10 minutes without use. (In the meantime the display is still displayed) otherwise you can turn it off by pressing the ON/OFF button. If the temperature is above 44°C, the display shows HC; if it is below 32°C, the temperature shows LC.

Comparative chart of normal temperatures recorded on various parts of the body:

Ear	between 35,8°C and 37,6°C
Mouth/forehead*	between 35,5°C and 37,8°C
Underarm	between 34,8°C and 37,4°C
Rectum	between 36,4°C and 38,2°C

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Dear customer,
We thank you for buying a Terraillon product.
We hope it proves entirely satisfactory and recommend you read this instruction leaflet carefully.

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