# DOWNLOAD WELLNESS COACH – HEALTH APP

- 1. Go to the App Store (iOS) or Google Play (Android) to download the Wellness Coach Health app.
- 2. Once you have installed the app, follow the instructions to create your account.

# SYNCHRONISE YOUR BODY COMPOSITION ANALYZER

- 1. Enable Bluetooth in your smartphone settings.\*
- 2. When the screen of your bathroom scale has switched off, keep the button on the back of the product pressed down for 3 seconds to launch synchronisation: the screen will show a rotating symbol indicating that synchronisation is in progress.
- 3. In the app:
  - A. Open the side menu.
  - B. Click on "Synchronise my devices".
  - C. Click on "Start synchronisation".
  - D. Once the app has detected your device, click on "Synchronise all devices".
  - E. Select a profile and enter your name, then click on Confirm.
- 4. Your body composition analyzer is synchronised with your smartphone. If it fails, or if the synchronisation process lasts for over a minute, "E1" will appear on the screen of your device.
- 5. Click on « Access my data » to return to the dashboard.

# START MEASURING

- 1. Place your body composition analyzer on a flat, hard surface.
- 2. Open the Wellness Coach Health app then check that your smartphone Bluetooth is ON\*. For the first weighing, select your profile number by pressing the arrows when your device is switched off, press SET to confirm your profile number.
- 3. Step onto your body composition analyzer.
- 4. Once the weight has stabilised, the scale measures your body composition.
- 5. The Bluetooth symbol appears to show that the connection to the app was successful. From now on, you will be able to consult the data transmitted by your device from your smartphone on the Wellness Coach Health app.

For all the next weighings, be sure that your smartphone Bluetooth is ON and that the Wellness Coach – Health app is openned before stepping onto your body composition analyzer. It will recognize automatically which profile you are thanks to your weight measurement.

### \*For Android users:

From Android version 6.0 onwards, Google will ask you to allow the Wellness Coach - Health app to access GPS data to enable your Bluetooth. You must provide this authorisation then enable the GPS function to be able to use the features of your scale which require a Bluetooth connection.

# FIRST USE AND SYNCHRONISATION OF YOUR ACTIVITY BAND

- 1. The first time you use your activity band, you must connect it to a USB port to activate it (see chapter "how to charge the device".
- 2. Open the Terraillon Wellness Coach app.
- 3. When selecting "Synchronise my devices → Activity band → Activi-T Smart", you have two options for synchronizing your activity band. Choose whichever method you prefer.

### a. OPTION 1: SCAN THE QR CODE PROVIDED

When the app asks you, scan the QR code that can be found in the box of your activity band (you can also find it on page 3 of the instruction manual provided in the box).

To avoid losing the QR code provided, make sure to stick it somewhere you will be able to find it easily. If you lose the QR code, you can still be able to synchronise your wristband using option 2.

#### **b. OPTION 2: FIND THE PRODUCT USING BLUETOOTH**

A 6-digit code will be sent to your wristband and displayed on the screen. Type this code into the app. Make sure to be quick, as the code will expire after one minute.

4. Once synchronisation is complete, you can access your dashboard by pressing the button "Access my data".

### HOW TO WEAR YOUR ACTIVITY TRACKER

- 1. Insert each part of the wristband onto the two ends of the pedometer.
- 2. Place the activity tracker around your wrist and adjust the size.
- 3. It is recommended to wear your wristband on your right wrist.

### HOW TO VIEW THE DATA

If the screen is off and there is enough battery, simply turn your wrist to turn it on. The screen will then display the following in this order: time and date, number of steps taken, calories burned, and distance covered. Then the screen will turn off.

Note: the data displayed on the screen will be reset at midnight. It is still saved on the device for 7 days if it has not been synchronised with the app.

When you receive a call on your smartphone, your wristband will alert you by vibrating and a notification will appear on the screen.

You can also set an alarm via the app in the "my sleep" section. At the alarm time, the alarm icon will appear and your wristband will vibrate.

### HOW TO CHARGE THE DEVICE

When the battery is almost or completely empty, insert the USB part of the device into a USB port (computer or other device with a USB port). Connect it the right way round (gold contacts) so that the battery icon starts flashing on the screen.

### WATER RESISTANCE

Your activity tracker conforms to the IP68 standard. This means that it is dust resistant and you can wear it in the shower or at the swimming pool. The wristband can withstand prolonged immersion up to a depth of 30 meters.

# DOWNLOAD WELLNESS COACH – HEALTH APPEN

- 1. Gå til App Store (iOS) eller Google Play (Android) for at downloade Wellness Coach Health appen.
- 2. Når du har installeret appen, følger du instruktionerne for at oprette din konto.

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