DOWNLOAD WELLNESS COACH – HEALTH APP

1. Go to the App Store (iOS) or Google Play (Android) to download the Wellness Coach - Health app.
2. Once you have installed the app, follow the instructions to create your account.

SYNCHRONISE YOUR BODY COMPOSITION ANALYZER

1. Enable Bluetooth in your smartphone settings.*
2. When the screen of your bathroom scale has switched off, keep the button on the back of the product pressed down for 3 seconds to launch synchronisation: the screen will show a rotating symbol indicating that synchronisation is in progress.
3. In the app:
 - A. Open the side menu.
 - B. Click on "Synchronise my devices".
 - C. Click on "Start synchronisation".
 - D. Once the app has detected your device, click on "Synchronise all devices".
 - E. Select a profile and enter your name, then click on Confirm.
4. Your body composition analyzer is synchronised with your smartphone. If it fails, or if the synchronisation process lasts for over a minute, "E1" will appear on the screen of your device.
5. Click on « Access my data » to return to the dashboard.

START MEASURING

1. Place your body composition analyzer on a flat, hard surface.
2. Open the Wellness Coach - Health app then check that your smartphone Bluetooth is ON*. For the first weighing, select your profile number by pressing the arrows when your device is switched off, press SET to confirm your profile number.
3. Step onto your body composition analyzer.
4. Once the weight has stabilised, the scale measures your body composition.
5. The Bluetooth symbol appears to show that the connection to the app was successful. From now on, you will be able to consult the data transmitted by your device from your smartphone on the Wellness Coach – Health app.

For all the next weighings, be sure that your smartphone Bluetooth is ON and that the Wellness Coach – Health app is opened before stepping onto your body composition analyzer. It will recognize automatically which profile you are thanks to your weight measurement.

**For Android users:*

From Android version 6.0 onwards, Google will ask you to allow the Wellness Coach - Health app to access GPS data to enable your Bluetooth. You must provide this authorisation then enable the GPS function to be able to use the features of your scale which require a Bluetooth connection.

DOWNLOAD WELLNESS COACH – HEALTH APPEN

1. Gå til App Store (iOS) eller Google Play (Android) for at downloade Wellness Coach – Health appen.
2. Når du har installeret appen, følger du instruktionerne for at oprette din konto.

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LAST NED APPEN WELLNESS COACH – HEALTH

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LADDA NED WELLNESS COACH – HEALTH-APPEN

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