Webcoach easyview

USER MANUAL | GUIDE D'UTILISATION

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iPad® 3/4/Air iPad® Mini

Android 4.3 Bluetooth Smart 4.0 / Bluetooth Smart Ready

If you can not find the app on Google Play from your Android smartphone, send us an email to service consommateurs@terraillon.fr with the model name of your smartphone.

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iPhone is registered trademark of Apple Inc. Android is a trademark of Google Inc. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. Other trademarks and trade names are those of their respective owners.

Webcoach easy-view

EN	TERRAILLON CONNECTED BODY COMPOSITION ANALYZER	4
FR	IMPÉDANCEMÈTRE CONNECTÉ TERRAILLON NOTICE D'UTILISATION	. 9
NL	TERRAILLON LICHAAMSVETWEEGSCHAAL HANDLEIDING	14
IT	BILANCIA MASSA GRASSA TERRAILLON MANUALE DI ISTRUZIONI	19
ES	BÁSCULA DE MASA GRASA DE TERRAILLON MANUAL DE INSTRUCCIONES	24
DE	TERRAILLON KÖRPERFETTANALYSEWAAGE BEDIENUNGSANLEITUNG	29
PT	BALANÇA DE MASSA GORDA TERRAILLON MANUAL DE INSTRUÇÕES	34

Thank you very much for selecting the Terraillon Glass Body Fat Analyzer WEB COACH EASY VIEW.

SAFETY INFORMATION

Don't use this product if you wear a pacemaker or implantable electronic device • Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening.

FIND MORE INFORMATION ON

http://www.terraillon.com

OVERVIEW

INITIAL START-UP

A Insert batteries in the scale (3xAAA) and in the screen (3xAAA)

B Press " \blacktriangle " button to change the numeral then press "SET" button to confirm. After confirming [DATE], [KG/LB/ST], [°C/°F] and [TIME] the LCD will display all

the information Please note that the time, date and your profile information (height, age, gender) will be updated automatically in the product when you synchronize it with the application in the following steps.

Press "MODE" button 3 sec to change the settings and use " $\blacktriangle \nabla$ " and "SET" to change them

Press "MODE" to use TIMER function(1/Countdown timer; 2/Toothbrush timer; 3/Shower timer). The screens «MEM?» and «- - - -» let you track your history. We recommend using your application to track your data.

D Download the Terraillon Wellness Coach App

E Create an account

F Turn on Bluetooth on your Smartphone (Setting > Bluetooth > ON)

G Hold 3sec "SET" button on the screen to start pair-up

It's important to bind your account with the scale for automatic recognition of user during weighing. You can pair the scale with 8 smartphones maximum. To do that, just repeat steps B to F on an other smartphone and select another profile number. To bind an other user on the same smartphone, please log in with the other user account and repeat steps E and F. Then select another profile number. Before you first weighing, please select your profile number on the scale.

If two profiles have a similar weight, the scale will ask you to confirm your profile number. On the scale, press the left arrow \blacktriangle if your profile number is the one on the left or press the right arrow $\mathbf{\nabla}$ if your profile number is the one on the right. On the below example, if you are profile 3, you have to press the left arrow \blacktriangle .

START MEASURING

Before your first weighing, it is necessary to choose a profile. Touch the SET Key, press up or down key to choose the profile Pn, and press SET to confirm. If you weigh yourself and your weight has changed more or less 2 kg compared to your previous weighing, please remember to select your profile before you weigh.

- A The scale will automatically switch on as you step on the platform barefooted. Stand still and keep full contact with the electrodes until the LED stops blinking ullu

Data ready to transmit to the App: If SUCCEED, the symbol disappears; If FAIL, the symbol Remains

The symbol Remains, successful Bluetooth Connection

WELLNESS COACH APP

Download the Terraillon Wellness Coach user manual to know more about our app

DELETE THE USER INFORMATION

A Delete single user information.

Touch SET key to enter the mode. Touch the down or up key to select the Person number (Pn with n = 1 to 8), touch the SET to confirm. When display the 0.0/ Pn (n is the profile number you want to delete), touch the up and down keys simultaneously more than 5 seconds. Then the product will delete person's information in the selected profile, and display "CLr" associated to the selected profile "Pn", and then back to the status of displaying the data, time and temperature.

B Delete All user information.

When display the 0.0, touch the up and down keys simultaneously more than 5 seconds. Then the product will delete all the person's information, and display "Clr A". Then back to the status of displaying the data, time and temperature.

MAINTENANCE

Use water only. DON'T use propellant or other chemicals.

SCREEN EASY VIEW

The product has Overload Instruction.

Low batteries on the weighing platform.

Low batteries on the LCD control terminal.

Body fat measurement error.

WARRANTY & ENVIRONMENT PROTECTION

This product is guaranteed for 2 years against material and manufacturing defects. During this period, such defects will be repaired free of charge (proof of purchase must be presented if the scale is under warranty). This warranty does not cover damage resulting from accidents, improper use or negligence. If you have a complaint please contact the store where you purchased the scale first.

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

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Terraillon SAS France & Headquarters
Terraillon UK Ltd
Terraillon Asia Pacific Ltd
Terraillon Corp USA contact@terraillon.fr