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
TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

Dear Customer,

Thank you for purchasing this product. We hope you will find it satisfactory. Please read these instructions carefully to take full advantage of your product.

NOTES ON SAFETY

Please read this section carefully to familiarize yourself with features and operations before using the unit.

- The warning signs and the sample icons shown here are listed in order for you to use this product safely and correctly as well as to prevent product damage, risk and injury to you or others.

 PRECAUTION WARNING	Indicate the right condition to use the product and prevent damage risk and injury.
 IMPORTANT FOR USERS	Indicate the important elements users should read before using the product.
 CARE AND MAINTENANCE	Indicate matters in which the possibility of damage may happen as a result of incorrect handling and improper maintenance.

PRECAUTION WARNING

- Do not use the scale on subjects who have a pacemaker or implantable electronic device. It may cause the device to malfunction. When in doubt, please consult your physician.
For pregnant women and people who have body implants such as artificial limbs, contraceptive devices, metal plates or screws, the results given by the Body Fat Analyser may be inaccurate. However there is no health risk associated with using the Body Fat Analyzer.
- When in doubt, please consult your physician.
- Do not disassemble the scale as incorrect handling may cause injury.

IMPORTANT FOR USERS

- Make sure use only the type of battery stated (see “Section Insert the Battery”).
- The «Athlete» mode is only applicable to those aged 18 or over.
- Body fat percentage estimates will vary with the amount of water in the body, and can be affected by dehydration or over-hydration due to such factors as alcohol consumption, menstruation, illness, intense exercise, etc.
- Do not use on pregnant women. The result is inaccurate.
- For body fat and body water estimates, subject must always be barefoot.


CARE AND MAINTENENCE

1. Do not disassemble the scale other than replacing the battery; it contains no user serviceable parts. Damage to the scale may occur as a result of improper handling.
2. Remove battery when the scale is not used for a prolonged period of time.
3. Clean the scale after use with a dampened cloth. Do not use solvents or immerse the unit in water.
4. Avoid excessive impact or vibration to the scale, such as dropping it onto the floor.


TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

PRODUCT OVERVIEW

- ① LCD Display
- ② Functions keys
- ③ Electrodes for body composition analysis
- ④ Unit button
- ⑤ Battery compartment

LCD OUTLOOK & FUNCTION KEYS

PREPARATION BEFORE USE

Insert the Batteries

Open the battery cover on the back of the scale unit. Insert the batteries and/or pull the plastic tape to initialize the batteries the first time. Be sure that the polarity of the battery is set properly for proper function (the polarity as indicated inside the battery compartment cover). If you do not intend to use this unit for a long period of time, it is advisable to remove the batteries before storage.

Installation & Initialisation

Place the scale on a hard, flat surface. The scale needs to be initialized at the first time of use, each time it has been moved or after battery replacement. Press briefly and sharply with your foot in the middle of the platform. The display shows "0.0kg". Your scale is now ready for use.

If your scale has a conversion switch, you can change the measuring units from metric to imperial or vice versa (i.e KG/ST/LB) by simply pressing the conversion switch located on the under-side of the scale.

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

INSTRUCTION MANUAL

Set up personal data(gender/age/stature) before using, it can set up 8 users' data.

1. Hold 3s "set" to enter the mode of setup, the flashing Pn (n can be 1~8) means the memory of the data.

Press ▲ & ▼ to change number, then press "SET" to confirm.

2. Set up the gender: press UP, DOWN (can select female, male), then press SET to confirm.

3. Set up the stature: Press UP, DOWN, and then press SET to confirm.

4. Set up age: Press UP, DOWN, and then press SET to confirm.

5. Set up the physical Activity Level (AL) from 1 to 5: Press UP, DOWN, and then press SET to confirm.

The following table is a reference for your Physical Activity Level (PAL) according to the World Health Organization based on human data:

Lifestyle	Activity Level (AL)
No physical activity - Chair-bound or bed-bound	AL 1
Seated work with no option of moving around and little or no strenuous leisure activity (short walks, light gardening...)	AL 2
Seated work, moderate leisure activity (30 minutes 2 to 4 times a week)	AL 3
Standing work, moderate to high physical efforts (30 minutes 4 to 6 times a week)	AL 4
Strenuous work or extremely high physical activity (at least 2 hours of intensive training 3 times a week) = Athlete mode	AL 5

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

6. After setting up of the data, when LCD displays 0.0, stand on the scale platform directly, the weight value will be shown after a few seconds.

7. After the weight is stable, the flashing “0” means that the scale is measuring the fat, the “0” signal moves from left to right for one times, and then shows the data.

Measuring process

Weight evolution since last measurement (kg)

Body Mass Index
(See page 19 for meaning)Body Fat (%)
(See page 19 for meaning)

Body Fat evolution since last measurement (points)

Body water (%)
(See page 20 for meaning)

Body Water evolution since last measurement (points)

Body muscle (%)
(See page 21 for meaning)

Body Muscle evolution since last measurement (points)

Daily Calorie requirement
(See page 22 for meaning)

NB: The evolutions are displayed only from the second weighing.

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

Last measurement review

1. When the scale is off, press UP or DOWN key to go to the desired memory number and press SET to confirm the user number to review last measurements.
2. Previous weight and previous Body Mass Values will be shown, then scale will go to off mode.

Automatically identify the user measurement

1. When the user stand on the scale directly, the scale will turn on and lock the weight.

2. When the weight is stable, display moving "0" means it is measuring:
 - If it can identify the user successfully, it will display the measurement.
 - If it detects many similar weight, it will display the relative user No.. Use UP key to choose the left user No. and use DOWN key to choose the right user No. After finished choosing the user, it will display the measurement and turn off.
 - If you don't choose the user, it will display the user No. several seconds and turn off.

Note: If 3 or 4 users have the same weight, screen is changing to display 2x2 users.

3. If it occurred error or it can't identify the user, it will display the present weight and turn off.


Delete the user information

When showing 0.0kg from the profile setting mode or last measurement review mode, touch UP and DOWN keys simultaneously for 5 seconds. The profile setting of the selected user number will be reset and its memory will be deleted. "CLr" will be shown for 3 seconds and go for 0.0 kg

WARNING INDICATORS

Overload: remove weight immediately Incorrect

Batteries failure:
Replace batteries

Impedance
measurement error

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

PRODUCT SPECIFICATION

- Bioelectrical Impedance Analysis (BIA) technology for body fat and body water, and muscle mass estimation (4 Poles solution) (BF: 0.1 %, BW: 0.1 %, BM: 0.1 kg, MM: 0.1 %)
- Min/Max capacity: 10/150kg
- BMI calculation
- DCR functionality
- Measurement difference
- 8 User Memories
- 3 touch keys operation (UP/SET/DOWN)
- Auto on and off function
- Overload and Low battery indication
- Age range from 10 to 85 years
- Physical activity levels from 1 to 5
- Height range from 75 to 225cm (2'5.5" to 7'4.5")
- User recognition with +/-2 kg ; +/-4,4 lbs
- Unit change (kg / st – lb / lb) by pushing the unit button behind the scale.

IMPORTANT INFORMATION TO KNOW BEFORE USING YOUR BODY FAT SCALE

1. What is BMI (Body Mass Index)?

- For adults from 19 to 99 years*:

It is the ratio between a person's weight and height. This index is used to assess a person's body fat. Calculation:
 $BMI = \text{weight (kg)} / [\text{height (m)}]^2$

* For adults over 65, the BMI classification can sometimes be inaccurate. However, using the Body Control does not in any way represent a health risk.

> 30	Between 25 and 29,9	Between 18,5 and 24,9	< 18,5
Obese	Overweight	Normal	Underweight

2. What is bio-electrical body mass analysis?

Bio-electrical body mass analysis is the resistance created by the human body when an electric current passes through it.

All body fat analysis bathroom scales work fundamentally on the same principle, which consists of measuring the bio-electrical impedance of the human body, and combining it with the subject's morphometric parameters (height, weight, gender, age, etc.) so as to calculate his/her pre-supposed amount of body fat, based on clinical measurements.

3. Why is it important to monitor percentage body fat (%BF)?

Weight is a useful indicator but not sufficient to monitor your health. Weight change in itself does not indicate whether it was the weight of body fat or muscle that had changed. In weight management, it is desirable that muscle mass be maintained while body fat is lost. Thus, monitoring the percentage fat of the body is an important step toward successful weight management and body health.

The optimal %BF of an individual varies according to age and gender.

The table as follows may be used as a guide:

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

Standard for Men

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	< 13 %	< 14 %	< 16 %	< 17 %	< 18 %
Optimal	14-20 %	15-21 %	17-23 %	18-24 %	19-25 %
Moderate	21-23 %	22-24 %	24-26 %	25-27 %	26-28 %
High	> 23 %	> 24 %	> 26 %	> 27 %	> 28 %

Standard for Women

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	< 19 %	< 20 %	< 21 %	< 22 %	< 23 %
Optimal	20-28 %	21-29 %	22-30 %	23-31 %	24-32 %
Moderate	29-31 %	30-32 %	31-33 %	32-33 %	33-35 %
High	> 31 %	> 32 %	> 33 %	> 34 %	> 35 %

4. How is percentage body fat (%BF) estimated?

%BF is measured by a method called Bioelectrical Impedance Analysis (BIA). With the BIA technology, a low intensity electrical signal is sent through the body. The signal is very low and causes no bodily harm. Depending on the amount of body fat of the individual, the electrical signal will flow with a different degree of difficulty. The difficulty with which a signal flows through the body is called electrical impedance. Hence, by measuring the electrical impedance and applying to the data a proprietary algorithm, %BF can be estimated.

Please be reminded that the %BF estimated with the scale represents only a good approximation of your actual body fat. There exist clinical methods of estimating body fat that can be ordered by your physician.

5. Why is it important to monitor percentage Total Body Water (%TBW) in the body?

Water is an essential component of the body and its level is one of the health-indicators. Water makes up approximately between 50-70% of the body's weight. It is present proportionally more in lean tissue compared to fat tissue. Water is a medium for biochemical reactions that regulate body functions. Water provides form to cells; helps to maintain body temperature; provides moisture to skin and mucosa; cushions vital organs; lubricates joints and is a component of many body fluids. The amount of water in the body fluctuates with the hydration level of the body and state of health. Monitoring the level of body water can be a useful tool for one's health maintenance. Similar to body fat estimation, the %TBW function provided in this scale is based on BIA.

The estimated %TBW may vary according to your hydration level, that is, how much water you have drunk or how much you have sweated immediately prior to the measurement. For better accuracy, avoid fluctuation in hydration level prior to the measurement. The accuracy of the scale in estimating TBW will also decrease with individuals suffering from diseases that tend to accumulate water in the body.

The optimal %TBW of an individual varies according to age and gender.

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

The table as follows may be used as a guide:

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Please be reminded that the % TBW estimated with the scale represents only a good approximation of your TBW. There exist clinical methods of estimating body fat that can be ordered by your physician.

6. When should I use the scale's body fat and total body water functions?

For maximum accuracy and repeatability, it is recommended that the scale's body fat and total body water functions be used at approximately the same time of the day, e.g. before breakfast in the morning. It is also a good practice to avoid swings in hydration level of the body prior to the measurement. Establishing your own baseline value of %BF and %TBW and track their changes is better than merely comparing your %BF and %TBW value to the population's normal value.

7. About muscle mass

Why is it important to monitor your muscle mass?

One of the key problem associated with getting older is a gradual decline in muscle mass. The scale gives you an estimate of your muscle mass with regards to your personal parameters (height, age, sex) as well as your total weight and body fat percentage.

There are no recommended results for muscle mass. However, they can help you monitor your body composition and its changes: as you increase your activity, your weight may not change but you may transform your body fat into muscles and change your body shape.

8. Why is the Athlete Mode (AL 5) necessary in a Body Fat Analyzer?

It has been found that body fat estimation using BIA could overestimate the percentage body fat of adult elite athletes. The physiological variation of athletes in bone density and level of hydration are two of the reasons said to account for the difference. The Athlete mode is selectable only for adults of 18 years of age or older.

9. Definition of an Athlete

The general consensus among researchers is that a quantitative dimension could be used in defining an athlete. For example, an athlete could be defined as a person who consistently trains a minimum of three times per week for two hours each time, in order to improve specific skills required in the performance of their specific sport and/or activity. Only in this case you can select Activity Level 5.

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

10. About Daily Calorie Requirement (DCR)

This is an approximate estimation of the daily calorie intake to maintain your present weight and lean body mass at your present physical activity level.

What are calories?

In medicine, the quantity of energy contained in food is measured in kilocalories (one kilocalorie = 1,000 calories). The human body uses this energy to live and breathe and go about its everyday activities. We can use the notion of kilocalorie (Cal) or kilojoule, which is now the international unit. Diets however are counted in calories. Just remember that 1 Cal = 4.18 kilojoules (kJoules, or kJ).

All food provides calories and all calories provide energy. The components of food (nutriments: protein, fat, carbohydrate) cover our energy needs.

Nutriments (for one gram)	Calories	Kilojoules
Carbohydrate	4 kilocalories	17 kilojoules
Protein	4 kilocalories	17 kilojoules
Fat	9 kilocalories	38 kilojoules

i.e. 1 g of proteins provides 17 kJ (4 kCal)

Other nutriments (Water, Vitamins, Minerals and Fibre) do not provide energy. Food and drink are composed of a combination of Protein, Fat and Carbohydrate. The energetic value of food therefore corresponds to the sum of calories provided by each nutriment. It is always given per 100 g of food.

Which factors influence the calculation of a person's energy requirements?

Several factors influence the calculation of a person's energy requirements. They correspond to the different origins of the body's energy expenditure.

- a- The basic metabolism: it is the minimum energy needed for everyday life. Even when resting, the human body is very busy: the heart beats, the brain sends its messages, the muscles work... the Basic Metabolism represents 60 to 70% of energy expenditure.
- b- Age: children and teenagers burn more energy per kilogram of body weight than adults. Energy requirements increase until around the age of 20 then they start to reduce by 2% per decade until the age of 60 and then by 10%.
- c- Gender: men spend more energy than women. Their body contains less body fat and 10 to 20% more muscle than women of the same age and height. Generally speaking, the energy requirements of men are 5 to 10% higher than those of women.
Exceptions: a pregnant or breastfeeding mother needs a higher energy intake for the good health of mother and baby.
- d- Thermoregulation: it corresponds to the energy needed to maintain the body at its physiological temperature of 37°C. A person living in a cold environment will burn more energy. Modern living conditions (central heating, warm clothes) reduce this consumption. This variation factor has become negligible.
- e- Thermogenesis: this is the energy burnt due to food intake, digestion of food and storage of nutriments. It represents 10% of total energy expenditure.
- f- Muscular work: an active person burns more calories than a sedentary person. Washing, dressing, gardening, walking or sport are all activities that consume the energy provided by our body. Expenditure depends on the type of activity and its intensity (Activity levels from 1 to 5).

TERRAILLON BODY FAT SCALE – INSTRUCTION MANUAL

Why is it important to monitor your Daily Calorie Intake?

Daily Energy Expenditure is defined by the OMS as “the quantity of energy needed to compensate for energy expenditure and ensure a size and body composition compatible with long-term good health and physical activity adapted to the economic and social context”. The most consequential variation factors are the Basic Metabolism and the Level of Physical Activity.

The balance between expenditure and energetic intake ensures the stability of weight and good health. Intake lower than expenditure leads to weight loss by the reduction of body fat and muscle and, in certain cases, nutritional deficiencies.

Inversely, intake that is higher than expenditure leads to weight gain by the storage of any excess as fat. Obesity and its pathological consequences can appear.

Nutritionists, dieticians and other healthcare professionals agree that a healthy diet and reducing fat intake are the two essential factors for healthy and durable weight loss. Losing weight by reducing the daily calorie intake helps develop knowledge of what is good for your body. Healthy and durable weight loss is possible by following a balanced diet and by reducing the daily intake of calories until it is lower than energy requirements (you burn more calories than you eat).

For example, for an adult aged 20 to 40, with the usual activities of the majority of the population, the recommended energy intake is 9200 kJ (2200 kcal) for women and 11300 kJ (2700 kcal) for men.

Regular weight watching and a very varied diet are two simple means of controlling balanced food intake. Be careful: dieting should be monitored by your doctor.

PROTECT THE ENVIRONMENT

Dispose of used batteries in a designed container so that they can be collected and recycled. Do not mix different types of batteries. Do not mix new and used batteries. After life, dispose of this product at a designated waste recycling point.

GUARANTEE

This product is guaranteed against defects in materials or manufacturing. During the guarantee period, any such defects will be replaced free of charge (proof of purchase must be shown in the case of a complaint under guarantee). The guarantee does not cover damage resulting from improper use or negligence. In the event of a claim, first contact the shop where you purchased the product.

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Terraillon SAS France & Headquarters
Terraillon UK Ltd
Terraillon Asia Pacific Ltd
Terraillon Corp USA