WEBCOACH SMART

USER MANUAL | GUIDE D'UTILISATION

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(🖾 🖏

iPhone is registered trademark of Apple Inc. Android is a trademark of Google Inc. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. Other trademarks and trade names are those of their respective owners.

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Any information provided by this device is in no way meant to treat, cure or prevent any disease or illness from happening.

FIND MORE INFORMATION ON

http://www.terraillon.com

PRODUCT OVERVIEW

FIRST USE

- Insert the batteries (3xAAA).
- You will need your smartphone to install your bathroom scales.
- Download the Wellness Coach Health app, then open it.

- Follow the instructions that appear on your smartphone. Ensure you enter all the information requested.
- Your smartphone must remain near to the bathroom scales throughout the installation process.
- If you encounter any problems, please consult the help centre:

START MEASURING

When the synchronization with your scale has been completed and the scale display is off, please select your profile number on the scale: use the arrows ▼▲ on your scale and press SET to confirm your profile. Step on the scale, the symbol \$ remains fixed on the display: the connection with the application succeeded.

Tip: After synchronizing with your scale, you can mount directly on your scale before the screen goes out. In this case, the weight data will be assigned to the selected profile during synchronization.

If the symbol does not appear, please ensure that your Wellness Coach- Health app is open, that Bluetooth is enabled on your smartphone and that you have selected the correct profile on the scale. If the symbol appears on the display of your scale before weighing, it means that the data from your previous weighing has not been transferred and has been saved in your scale. If the symbol a disappears from the display of your scale after weighing, it means that all the data has been transferred.

If two profiles have a similar weight, the scale will ask you to confirm your profile number. On the scale, press the right arrow $\mathbf{\nabla}$ if your profile number is the one on the right. In the example below, if you are profile 1, you must press the left arrow $\mathbf{\Delta}$.

Be sure to always weigh yourself on a hard and flat surface.

MAINTENANCE

Use water only. DON'T use propellant or other chemicals.

WARRANTY & ENVIRONMENT PROTECTION

This product is guaranteed for 2 years against material and manufacturing defects. During this period, such defects will be repaired free of charge (proof of purchase must be presented if the scale is under warranty). This warranty does not cover damage resulting from accidents, improper use or negligence. If you have a complaint please contact the store where you purchased the scale first.

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice.

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Terraillon SAS France & Headquarters
Terraillon UK Limited
Terraillon Asia Pacific Ltd
Terraillon Corp USA